

Product Manual

LentiBOOST[®] Transduction Enhancer

Shipped at room temperature Store at -20 °C

FOR RESEARCH USE ONLY (LentiBOOST®-P)

Revision 01/2022

www.sirion-biotech.com



1. **PRODUCT DESCRIPTION**

LentiBOOST® enhances the uptake of lentiviral vectors into mammalian cells. Compared to commonly used transduction enhancers like Polybrene or Protamine sulfate it does not negatively affect cell viability or cell growth.

Furthermore, for human CD34⁺ hematopoietic stem cells and progenitor cells LentiBOOST® was shown to maintain the "stem-cell like" potential of the cells after lentiviral transduction (Hauber et al. 2018). Moreover, it shows high transduction rates while limiting the vector copy number (VCN < 5) that reaches each individual cell even at low MOI (< 10).

For human T-cells it was shown that the use of LentiBOOST® for lentiviral transduction results in a 2-fold increase in transduction rate, has no cytotoxic effect on the cells and does not affect cytokine production after antigen encounter (Simon et al. 2019).

Additionally, the use of LentiBOOST® can increase vector copy number (VCN) up to 4-fold (Hauber et al. 2018).

LentiBOOST® is therefore an excellent tool for enhancing lentiviral transduction of sensitive primary cells including hematopoietic stem cells and T-cells.



Fig. 1: Lentiviral transduction efficiencies in H1299 48h after transduction with lentivirus LV-CMV-GFP. Protocols according to the instructions of this manual



2. Material

LentiBOOST® 100 mg/ml aqueous solution

3. Storage

Store at -20°C.

4. Thawing and handling

It is recommended to thaw LentiBOOST® at room temperature <+25°C.

Before opening the tube spin down briefly to remove any liquid from the lid. For further use aliquot LentiBOOST®-P using aseptic technique and store at -20°C. LentiBOOST®-P can be stored at +4°C for 1 month. No loss of functionality* was observed after 5 freeze-thaw cycles.

*based on an internal, standardized assay, no conclusion can be drawn to any customer specific cell types

PLEASE NOTE: the protocols given below are suggestions based on (published) customer data. Protocols have to be adapted depending on customer specific conditions.

5. Transduction Protocol for CD34+ Transduction

Preparation of human CD34+ HSC

CD34+ hematopoietic stem cells (HSC) are isolated from PBMC using standard protocols

Day 1: Seeding cells

CD34+ HSC are cultured according to standard protocols. Cells are incubated at 37°C in a humidified incubator containing 5% CO2.

Day 2: Transduction

The following day, CD34+ HSCs are counted and $1x10^6$ HSCs transduced in the presence of LentiBOOST® and Lentivirus particles in 24 well plates using MOI 10 (see note below). For an initial experiment it is recommended to use MOIs between 2-30 for transduction and to add LentiBOOST® at the standard concentration of 1 mg/ml (1:100) of the total volume (medium + virus). In a second experiment it is recommended to titrate LentiBOOST® in the range of 5 mg/ml – 0.1 mg/ml (1:20-1:1000) to determine the minimal active concentration.



Protocol steps:

• Calculate the volume of lentivirus needed (see Table 1)

					Amount of infectious lentiviral particles					
	Plate	CD34 Cell number/well	Total Volume (Medium+Virus)	Volume LentiBOOST®	MOI 30	MOI 15	MOI 10	MOI 5	MOI 2	MOI 0
2	24 well	1.00E+06	1 ml	10 µl (1:100)	3E+06	1.50E+06	1.0E+06	5.0E+05	2.0E+05	0

Table 1: Recommended volumes of medium and LentiBOOST®-P to be used for lentiviral transductions in an initial experiment. It is recommended to titrate LentiBOOST®-P in a second experiment in order to determine the minimal active concentration. Therefore, we recommend to dilute LentiBOOST® in the range of 1:20-1:1000 (5 mg/ml – 0.1 mg/ml). For other multiwell plates parameters have to be adjusted accordingly.

- Thaw lentivirus at +4°C
- Add the appropriate amount of LentiBOOST® according to Table 1 directly to the cells seeded the day before
- Add the amount of lentivirus according to Table 1 directly to the cells and mix carefully
- Incubate cells over night at normal cell culture conditions

Day 3: Medium exchange

- Aspirate medium from transduced cells and add appropriate amount of normal growth medium
- Note: The MOI was determined according to the following method: VSV-G pseudotyped lentiviral vector particles purified and concentrated by ultracentrifugation were resuspended in StemMACS™ HSC Expansion Media (Miltenyi Biotec), aliquoted and stored at -80°C.
- The titer of pseudotyped lentiviral particles was determined as fluorescence forming units per ml (ffu/ml). This required transduction of 293T cells grown in standard medium (DMEM 10 % FBS) in 24-well plates with different volumes of viral vectors in presence of 0.5 µg/ml protamine sulfate. The cells were then spinoculated at 600 x g for 10 min at room temperature. Medium was changed 8 h after transduction. At 72 h post transduction, cells were harvested and analyzed by flow cytometry for GFP expression. Samples that contained 5 to 25% GFP positive cells were used to calculate viral titers.

6. Transduction Protocol T-cells

Seeding and pre-stimulation of cells

• Primary T Cells are stimulated with anti-CD3 antibody, anti-CD28 antibody and interleukin-2 (IL-2) according to standard protocols. This step should be adapted according to individual needs.

Day 1: Transduction

- Thaw lentivirus at +4°C
- Prepare 500µl of culture medium and add LentiBOOST® at 1 mg/ml (1:100) as starting point. Concentrations of LentiBOOST® can be tested between 0.1 mg/ml and 5 mg/ml.
- Add viral vector to medium at desired MOI and mix gently. For an initial experiment it is recommended to use MOIs between 2-30.
- Pellet 10⁶ cells by centrifugation (cell number can be adapted according to needs)



- Mix the cell pellet with the prepared medium containing LentiBOOST® and viral vector
- Seed the cells into a 24 well plate
- Spinoculation (optional): Centrifuge cell culture plate for 90 min at 800 g at RT.
- Incubate the overnight at 37°C and 5% CO₂

Day 2: Medium change

• Exchange medium according to standard protocols

- Hauber, I., Beschorner, N., Schrödel, S., Chemnitz, J., Kröger, N., Hauber, J., & Thirion, C. (2018). Improving Lentiviral Transduction of CD34+ Hematopoietic Stem and Progenitor Cells. *Human Gene Therapy Methods*, hgtb.2017.085. https://doi.org/10.1089/hgtb.2017.085
- Simon, B., Harrer, D. C., Thirion, C., Schuler-Thurner, B., Schuler, G., & Uslu, U. (2019). Enhancing lentiviral transduction to generate melanoma-specific human T cells for cancer immunotherapy. *Journal of Immunological Methods*, 472(April), 55–64. https://doi.org/10.1016/j.jim.2019.06.015

Trouble	Possible reason	Solution			
	MOI used was too low	Use higher amounts of lentivirus up to MOI 50			
Low transduction		Include spinoculation step 800 g for up to 90 min at room temperature (in cell culture plates)			
efficiency	Cells are very hard to transduce	Increase concentration of LentiBOOST-P up to 5 mg/ml (1:20)			
		add Protamine sulfate at 5 µg/ml additionally to LentiBOOST(R)			
	Cells are sensitive to LentiBOOST®	Decrease concentration of LentiBOOST® to			
		e.g. 0.2 mg/ml, 0.1 mg/ml (1:500, 1:1000)			
		Try protocol without spinoculation			
l ow viability	Cells are sensitive to spinoculation	Reduce duration			
Low viability		Reduce velocity			
	Cells are sensitive to lentiviral vectors	Change medium 4 h after transduction or directly after centrifugation			

8. Troubleshooting:



Limited Use Label License Agreement

IMPORTANT INSTRUCTIONS – READ CAREFULLY: This Label License Agreement (this "Agreement") is the legal agreement between you (hereinafter "Licensee") and SIRION Biotech GmbH ("SIRION") (individually, a "Party" or together, the "Parties"), for poloxamers or products containing poloxamers ("Product") and use thereof, whether such Product is provided by SIRION, its affiliates, (sub)licensees, resellers, or distributors. If the Product is subject to multiple limited use label licenses, the terms of the most restrictive limited use label license shall control. Licensee's use of the Product shall constitute acknowledgment and acceptance of the terms and conditions of this Agreement. If Licensee is not willing to use the Product pursuant to the terms and conditions of this Agreement, please contact SIRION Biotech GmbH to return the unused and unopened Product for a full credit. Address: Am Haag 6, 82166 Gräfelfing, Germany Email: Licensing@sirion-biotech.com, Tel.: +4989700961999

1. Use of Products by Licensee. SIRION holds a patent (published as WO2013127964) for this substance class and licenses the technology. SIRION grants to Licensee a non-exclusive, non-transferable, limited right, subject to certain restrictions as set forth below, to use only the purchased amount of the Product for Research Purposes Only; provided that no license, right or permission is granted hereunder for a non-academic, for-profit, or commercial Licensee to use the Product in the Excluded Field. "Research Purposes Only" means internal basic research as conducted by Licensee for the sole benefit of Licensee. "Research Purpose Only" specifically excludes Screening, Contract Research or Contract Screening. "Excluded Field" means any use of the Product for *ex vivo* gene therapy for hemoglobinopathies. The Product and use thereof are protected under patent applications/patents (WO2013127964) and other intellectual property rights owned or controlled by SIRION or its affiliates. SIRION and its affiliates retain all rights not expressly granted herein, and there are no implied licenses or rights granted herein by estoppel, implication or otherwise.

2. Excluded Uses. The Product is <u>not</u> licensed, and Licensee may <u>not</u> use, the Product for any Commercial Purpose, Development Purpose or in the Excluded Field, and no licenses are granted expressly, by implication, by estoppel, or otherwise for any Commercial Purpose, Development Purpose or in the Excluded Field. "Commercial Purpose" means any commercial activity, including but not limited to: (a) any use, directly or indirectly, in manufacturing and/or production, or quality control; (b) any use to provide a service, information, or data; (c) any use for therapeutic, diagnostic, vaccine or prophylactic purposes; (d) any sale, resale, leasing, distribution, transfer, licensing, and/or importation of, whether or not for Research Use Only; and (e) any applications which require regulatory approval, including any *in vitro* diagnostic or therapeutic applications in humans, any *ex vivo* use to enable gene therapy, or any *in vivo* use in humans for any purpose. "Development Purpose" means any clinical activity following IND-enabling preclinical toxicological studies or equivalents thereof. Licensee acknowledges that Licensee shall not have the right to authorize any third party to use or sell the Product or derivatives of the Product which contain poloxamers protected by patent applications/patents and other intellectual property rights owned or controlled by SIRION or its affiliates. The foregoing amounts to a restricted license under applicable patent and other intellectual property rights.

3. Misuse of the Products and Indemnification terms. To the extent provided by law, Licensee will defend, indemnify and hold harmless SIRION and its affiliates, their managers, directors, officers, employees, sponsors, and agents (collectively the "Indemnified Parties") against any and all liability, loss, damage, claim or expense, including attorney's fees (collectively the "Indemnified Losses") arising out of or in connection with this Agreement, including, without limitation Indemnified Losses resulting from any use by the Licensee, Licensee's employees or students of the Product and any derivatives of the Product. Licensee will indemnify and hold harmless the Indemnified Parties against any and all Indemnified Losses resulting from, arising out of or relating



to: (i) product liability claims of any nature; (ii) claims arising from Licensee's failure to comply with all governmental regulations relating in any way to use or storage of the Product; (iii) Licensee's breach of this Agreement; and (iv) claims by a third party that Licensee, Licensee's representatives, personnel or students use of the Product infringes or violates any patent, copyright, trademark or property rights of such third party. This Agreement shall be interpreted and enforced in accordance with the laws of the State of New York, without regard to its conflict of law rules.

4. DISCLAIMER OF WARRANTY. NEITHER SIRION NOR ANY OF ITS AFFILIATES MAKES ANY WARRANTY, EXPRESS, OR IMPLIED WITH RESPECT TO THE PRODUCT, INCLUDING ANY WARRANTY OF MERCHANTIBILITY OR FITNESS FOR ANY PARTICULAR PURPOSE OR THAT THE PRODUCT DOES NOT INFRINGE ANY PATENT, COPYRIGHT, TRADEMARK, OR OTHER RIGHTS OR ANY OTHER EXPRESS OR IMPLIED WARRANTIES.

5. Term. This Agreement shall commence upon the receipt of the Product and, subject to Section 6, continue as long as the Product is in the possession of the Licensee.

6. Termination. Either Party may terminate this Agreement for any reason upon thirty (30) days' written notice to the other Party. The rights and obligations under Sections 3, 4, 6, and 7 shall survive any termination, expiration, or completion of this Agreement with respect to information generated and activities and events occurring prior thereto. Upon expiration or any termination of this Agreement, Licensee shall promptly destroy all remaining Product and all material made using the Product.

7. Entire Agreement and Assignability. This Agreement sets forth the complete and entire agreement of the Parties with respect to the Product and supersedes and terminates all prior agreements and understandings between the Parties with respect to the Product. No subsequent amendment or addition to this Agreement shall be binding upon the Parties unless reduced to writing and signed by the respective authorized representatives of the Parties. This Agreement shall not be assigned or otherwise transferred by Licensee.

Additional Rights

Licensee will promptly inform SIRION in writing at the address specified below in the event of any uses in conflict with those permitted by this Agreement. For information on obtaining rights for uses other than those permitted by this Agreement, please contact SIRION:

Business Development Licensing SIRION Biotech GmbH Street: Am Haag 6 Postal code: 82166 City: Gräfelfing Country: Germany Email: Licensing@sirion-biotech.com Tel.: +49 89 70096199-9

SIRION reserves the right to grant, or not to grant, permission for any such other uses.